

Protocol for chemical inducible expression of MBP & sample preparation

Materials:

- HgCl₂ store solution with the concentration of $10^{-3}M$
- Overnight bacterial culture or bacterial colonies;
- fuming nitric acid
- IPTG

1. Add 1ml of the overnight bacterial culture or pick a colony to 100ml of LB antibiotic medium, Incubate at 37 degree in a shaker till the OD₆₀₀ value reaches 0.6.
2. Add 1ml IPTG to the 100ml bacterial to induce the expression of Mental Binding Peptide (MBP) or Lpp-Ompa-MBP or Dsba-MBP.
3. Incubate the culture at 30 degree for 0.5 hour.
- 4 Supply HgCl₂ to the culture to different concentrations of $10^{-6}M$, $10^{-7}M$, $10^{-8}M$ and control.
5. Incubate the culture at 30 degree for 40 hours.
- 6 Harvest cells by spinning at 5000rpm for 5 min. Discard the supernatant, Wash with 6ml of ddH₂O for two times.
- 7 Dehydration the bacteria with freeze dryer, followed by measurement of the weight of cell pellet.
- 8 Resuspend the pelleted cells in 8ml fuming nitric acid, followed by microwave digestion. Resuspend sample to 25mL. Then sample preparation

is completed.