Cell Survivability Testing

(Prepared by Wisconsin 2010 iGEM team)

Parts used in this experiment			
Part Number	Function	Expression Type	Zip File
<u>BBa_K318500</u>	Produces Trascription Factor RcsB	Inducible - IPTG	500
BBa_K318501	Produces Trascription Factor RcsA	Inducible - IPTG	501
BBa_K318502	Produces Trascription Factor RcsA & RcsB	Inducible - IPTG	502
BBa_K200021	Empty Vector/Contol	Inducible - IPTG	NA

Materials

1M HCl

1M NaOH

LB

1M IPTG

Protocol for Acid Survivability Test

- 1. Grow cells overnight in LB.
- Inoculate a new sets of samples which have same starting OD (OD=0.05) by using the overnight culture. Add 1mM IPTG to induce.
- Harvest cells at OD between 1 to 2 Wait until the slowest one is at OD=1.5
- 4. Make LB+ITPG cultures that are OD=1.5, and the volume should between 1ml to 5ml.
- 5. Spin down the cells at 3000g for 2-3mins.
- 6. Re-suspend the cells in pH LB medium (pH=2, 4, 7)
- 7. Plate cell samples at time zero.
- 8. Place all the samples into the 37C shaker and let them grow for 4hs.
- 9. Remove the samples from the shaker, and record the OD.
- 10. Make a series dilution for each sample. 10E-3, 10E-4, 10E-5, 10E-6, 10E-7, 10E-8, 10E-9
- 11. Plate 20ul of each dilution of each sample on a CM mini-Petridish plate.
- 12. Put all the plates in to 37C incubator and let the cells grow overnight.



Figure 1 (Created by Wisconsin-Madison 2010 iGEM team)



Figure 2 (Created by Wisconsin-Madison 2010 iGEM team)



Figure 3 (Created by Wisconsin-Madison 2010 iGEM team)