

Protocol 1-6: Preparation of Competent Cell for Electro Transformation

General Protocol

- 1) Add 20 μ l of the overnight bacterial culture or pick a colony to 1 ml of LB antibiotic liquid medium, incubate at 37 degree in a shaker till the OD600 value reaches 0.4-0.6.
- 2) Put the tubes on ice to incubate for 5 min.
- 3) Pellet bacterial cells by 4 min centrifugation at 4000 rpm, discard the supernatant
- 4) Resuspend cells in 600 μ l of ice-chilled 0.1 M Calcium Chloride solution. Incubate on ice for 30 min.
- 5) Centrifuge for 4 min at 4000 rpm in a microcentrifuge tube, discard the supernatant.
- 6) Resuspend the pelleted cells in 100 μ l of ice-chilled 0.1 M Calcium Chloride solution. Incubate on ice.
- 7) Add 50 μ l of the prepared cells to each tube containing DNA sample, mix and incubate on ice for 30 min.
- 8) Transform subsequently as the transformation protocol.
- 9)

Tips

- 1) Make sure the cells are not left in the centrifuge at ambient temperature for more than 5 min as this will significantly decrease the transformation efficiency.
- 2) The rpm at centrifugation is not higher than 4000, as a high rpm may cause the lysis of cells.
- 3) Competent cells prepared with this protocol are suitable for direct use only. Freezing down and storage at -70°C is not recommended.
- 4) The culture can be kept at 4 degree for one week and used for preparation of competent cells, but culture stored longer than 10 ten days is not suitable for competent cells.

Reference

- 1) Sambrook J, Maniatis T, Fritsch EF. Molecular Cloning: a Laboratory Manual. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY, 3rd ed., 2001.