

# LB medium and LB agar gel

## Medium for cultivation of *E. coli*

### ◆ LB medium (1 L)

- 1; Add about 900 mL of distilled water to beaker.
- 2; Add reagents (Table) and stir.
- 3; Add distilled water up to 1 L and take LB medium to media bottle.
- 4; Autoclave for 20 min at 120°C

Table Composition of LB medium

Reagent	Volume
Bact Tryptone	10 g
Yeast Extract	5 g
NaCl	5 g
For 1 L LB medium	

### ◆ LB agar gel (1 L)

- 1; Prepare LB medium without autoclave (Steps 1-3 of 1L scale of LB medium).
- 2; Add 15 g of agar and stirrer bar.
- 3; Autoclave for 20 minutes at 120°C.
- 4; Stir and cool LB medium with agar, add appropriate antibiotic (50 µL).
- 5; Pour LB medium (Step 4) in plate and cool down in clean bench.

## Transformation

### Preparation of *E. coli* contain particular plasmid

- 1; Incubate frozen competent cell (DH5α) on the ice for a few minutes.
- 2; Add 1~5 µL of plasmid to competent cell (DH5α) on the ice.
- 3; Incubate for 20 – 30 minutes on the ice.
- 4; Incubate for 45 seconds at 42°C.
- 5; Add 1 mL LB medium and cultivate for 1 hour at 37°C.
- 6; Spread culture medium on LB agar plate with appropriate antibiotic.

## Plasmid extraction

### Preparation of plasmid extracted from *E. coli*

Solution	Solution composition	Attention
Solution 1	50 mM Tris-HCl (pH8.0), 5 mM EDTA, 20 mg/ml RNase A	Stare at 4°C
Solution 2	0.2 M NaOH, 1% SDS	
Solution 3	3 M Potassium acetate (pH 5.5 adjusted by acetate)	
Solution 4	4 M NaClO <sub>4</sub> (pH 6.0 adjusted by Tris )	Stare at 4°C
Bind mix	5 M GuSCN, 4%(w/v) Triton X-100, Silica (0.25 g/100 mL)	Shading

- 1; Pick up single colony from agar plate and cultivate it in 1.5 mL LB medium containing appropriate antibiotic (50 µg/mL) overnight at 37°C.
- 2; Move the culture medium to 1.5 mL tube.

- 3; Centrifuge for 5 seconds at 15,000×g and 4 °C and discard supernatant.
- 4; Add 100 µL Solution 1 to the pellet and resuspend and incubate for 3 minutes.
- 5; Add 100 µL Solution 2, invert tube gently 5 times and incubate for 3 minutes.
- 6; Add 100 µL Solution 3, invert tube 5 times and incubate for 3 minutes.
- 8; Add 200 µL Solution 4 and invert tube 5 times and centrifuge for 3 minutes at 15,000×g and 4 °C.
- 10; Take supernatant to new 1.5 mL tube and centrifuge for 3 minutes at 15,000×g and 4 °C.
- 11; Vortex Bind mix for 1 min and add 800 µL Bind mix to new 1.5 mL tube.
- 12; Add 400 µL supernatant after centrifugation (Step 10) to tube containing Bind mix (Step 11) and mix.
- 13; Incubate for 3 minutes, centrifuge for 3 seconds at 5,000×g and 4 °C, and discard supernatant.
- 14; Add 1 mL of 50% ethanol and resuspend.
- 15; Centrifuge for 3 seconds at 5,000×g and 4 °C and discard supernatant.
- 16; Repeat wash (Steps 14-15).
- 17; Dry pellet for a few minuet under a vacuum to remove residual ethanol.
- 18; Add 50 µL nuclease-free water or TE buffer and incubate for 3 minutes at 65°C.
- 19; Centrifuge for 3 minutes at 15,000×g and 4 °C.
- 19; 40 µL of supernatant into new 500 µL tube .

## Restriction enzyme digestion of DNA

### Cleavage of insert DNA from plasmid

- 1; Mix DNA and restriction enzyme (Table).
- 2; Incubate for 2 hours at 37°C.
- 3; Incubate for 10 minutes at 65°C.
- 4; Confirm the band of DNA by agar gel electrophoresis.

Table Composition of LB medium

Reagent	Volume
DNA plasmid	5 µL
10 × buffer	1 µL
Nuclease-free water	3 µL
<i>EcoRI</i> or <i>XbaI</i>	0.5 µL
<i>SpeI</i> or <i>PstI</i>	0.5 µL
Total	10 µL

## Agar gel electrophoresis

### Confirmation and separation of digested DNA

#### ◆ Preparation of agar gel

- 1; Add 1 g of agar to 100 mL of 1 × TAE.
- 2; Boil and stir until solution is clear.
- 3; Cool down, pour to gel form and set gel corm.
- 4; Incubate until gel dry out.
- 5; Stare gel in 1 × TAE.

Reagent	Volume	Composition
50 × TAE	20 mL	2 M Tris, 2 M Acetate, 0.5 M EDTA (pH 8.0)
Distilled water	980 mL	
1 × TAE	1,000 mL	

## ◆ Agar gel electrophoresis

- 1; Place agar gel and pour 1 × TAE in electrophoresis chamber.
- 2; Load DNA ladder and DNA sample mixed with loading dye on agar gel.
- 3; Electrophorese for 20 minutes at 100 V.
- 4; Stain gel by Sybr Safe™ Gel Stain (Invitrogen).
- 5; Visualize the band of DNA using UV light.
- 6; Confirm the length of digested DNA

## Gel purification

### Purification of DNA from agar gel

#### GENECLEAN® II Kit(Nal, glass milk, NEW Wash)/Qbiogene

- 1; Cut the objective band in the agar gel after electrophoresis and stain with SYBR Safe.
- 2; Put the gel including objective DNA into 1.5 mL tube and measure the mass of that.
- 3; Add 2.5-3 fold volume Nal solution into the tube (Gel: Nal =1 mg : 1 µL ).
- 4; Incubate the gel at 50°C for 5 minute.
- 5; Add 10 µL of glass milk and vortex.
- 6; Incubate for 5 minutes and vortex per a minute.
- 7; Centrifuge for 5 seconds at 15,000×g and 4°C and discard the supernatant.
- 8; Add the 500µL of New Wash and resuspend.
- 9; Centrifuge for 5 seconds at 15,000×g and 4°C.
- 10; Repeat wash (Steps 8-9).
- 11; Dry the pellet for 5-10 minutes under vacuum.
12. Add 20 µL of nuclease-free water and resuspend.
13. Centrifuge for 5 seconds at 15,000×g and 25°C.
14. Transfer supernatant including objective DNA into new tube.

## Ligation

### Ligation inset DNA and vector

#### DNA Ligation kit Ver 2.1(Solution I )/Takara

- 1; Mix the insert DNA, vector and solution I (Table).
- 2; Incubation at 16°C for 30 minute.
- 3; Transform *E. coli* with ligation sample.

Reagent	Volume
Solution 1	7 µL
Backbone plasmid	1 µL
Upstream insert DNA	4 µL
Downstream insert DNA	4 µL
Total	14 µL

## Colony PCR

Confirm of insert DNA in plasmid, directly using *E. coli* at PCR

**Table 1 Composition of reagent solution for colony PCR**

Reagent	Volume	Final conc.
2 × Go Taq® Green Master Mix (Promega)	5 µL	1 ×
Forward primer (100 µM)	0.1 µL	1 µM
Reverse primer (100 µM)	0.1 µL	1 µM
Nuclease-free water	4.8 µL	
Total	10 µL	

**Table 2 PCR program**

Temp. (°C)	Time (min)	Cycle
95	2	1
95	0.5	} 30
55	0.5	
72	2*	
72	5	1

\*: Adjust time to 1 min/kbp

- 1; Add 10 µL of reagent solution (Table 1) to PCR tube.
- 2; Pick up single colony from agar plate with tooth pick and sting replica plate (new LB agar plate).
- 3; Put and stir toothpick to reagent solution (Step 1).
- 4; Amplify insert DNA with PCR program (Table 2).
- 5; Electrophorese PCR sample with agar gel
- 6; Check the band and length of insert DNA and decide the colony with insert DNA

## Sequence analysis

### Identification of insert DNA

#### ◆ Preparation of PCR product

**Big Dye® Terminator Cycle Sequencing Kit Ver. 3.1 (Premix, Buffer) / Applied Biosystems**

- 1; Add reagent solution (Table 1) to PCR tube and amplify insert DNA with PCR program (Table2).

**Table 1 Composition of reagent solution for colony PCR**

Reagent	Volume	Final conc.
Premix	1 µL	
5 × Buffer	4 µL	1 ×
Primer (1 µM)	4 µL	200 nM
Plasmid	1 µL	1 µM
Nuclease-free water	10 µL	
Total	20 µL	

**Table 2 PCR program**

Temp. (°C)	Time (min)	Cycle
95	5	1
96	0.5	} 25
50	0.25	
60	4	
4	∞	

#### ◆ Purification of PCR product and sequence analysis

**Agencourt CleanSEQ® and 96 R ring Super Magnetic Plate® / Beckman Coulter**

- 1; Add 10 µL of Agencourt CleanSEQ® 10 µL to PCR product.
- 2; Add 62 µL of 85% ethanol, mix and incubate for 3 minutes.
- 3; Incubate for 3 minutes on 96 R ring Super Magnetic Plate® and discard supernatant.

- 4; Add 100  $\mu\text{L}$  of 85% ethanol and mix.
- 5; Incubate for 3 minutes on 96 R ring Super Magnetic Plate<sup>®</sup> and discard supernatant.
- 6; Repeat wash (Steps 4-5).
- 7; Dry for 10 minutes.
- 8; Add 40  $\mu\text{L}$  nuclease-free water and mix.
- 9; Transfer 30 $\mu\text{L}$  of clear sample into a new plate for loading on the detector.
- 10; Load sample on sequencer and analyze.