

## Protocol 1-8: Electro Transformation of Recombinant DNA

### General Protocol

- 1) Thaw electrocompetent cells on ice, keep on ice once thawed
- 2) Mix DNA( < 2 ul with low salt content to prevent arcing) with 40 ul cells in an eppendorf tube, keep on ice or alternatively, or add DNA to the electro-cuvette (pre-chilled) and then add the electrocompetent cells making sure that the DNA gets mixed well with the cells.
- 3) Transfer the cell/DNA mix into an electroporation cuvette
- 4) Set the pulser properly( See Tip 1)
- 5) Pulse the cells once; the voltage display blinks, and the gene pulser beeps
- 6) Quickly transfer 37 C SOC to cuvette, mix by gently pipetting up and down, and transfer SOC/cells back to culture tube; replace in 37 C H<sub>2</sub>O bath
- 7) Incubate 30 min. @ 37 C, followed by 30 min. @ 37 C/shaking
- 8) Plate cells on appropriate antibiotic

### Tips

- 1) the gene pulser should be set properly
  - time constant = 4.5 - 5.0 ms
  - resistance = 200 W
  - capacitance = 25 mFDfor 0.1 cm gap cuvettes, set the volts to 1.7 kV  
for 0.2 cm gap cuvettes, set the volts to 2.5 kV

### Reference

- 1) Sambrook J, Maniatis T, Fritsch EF. Molecular Cloning: a Laboratory Manual. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY, 3rd ed., 2001.
- 2) Robert F. Weaver. Molecular Biology, McGrawHill, 4th edition, 2007