

Colony PCR flu operon
8/4/2010

Prepare DNA Template

1. Take a stab of the frozen stock of *E. coli* K12 from the IGEM -80C frozen stock (not too much because too much template will lead to unspecific replication) and resuspend it in 50 uL of DI water
2. Mix primers to a 50 mM concentration
 - a. Multiply the nmol concentration from the IDT printout by 20
 - b. This is the volume of resuspension in ultrapure water in uL
 - i. Forward flu primer: $24.5 \text{ nmol} \times 20 = 490 \text{ uL}$ ultra pure water
 - ii. Reverse flu primer $46.9 \text{ nmol} \times 20 = 938 \text{ uL}$ ultra pure water
 - c. Done previously
3. Turn on PCR machine to instant incubate at 98 C for the first initial denaturation step
4. Place PCR tubes for run on ice
5. PCR reaction KEEP EVERYTHING ON ICE
 - a. Vortex all tubes before starting to make sure everything is well mixed
 - b. For a single reaction mixture (20 uL)
 - i. 14.5 uL Ultra pure water
 - ii. 4 uL of 5x phusion master mix
 - iii. 0.4 uL 10 mM dNTP
 - iv. 0.25 uL of primer A (0.5 uM concentration from 40 uM stock)
 - v. 0.25 uL of primer B (0.5 uM concentration from 40 uM stock)
 - vi. 0.4 uL of DNA template (from step 1)
 - vii. 0.2 uL Phusion DNA polymerase
 - c. Combine all ingredients but polymerase for 5 samples in master mix (100 ul)
 - i. 72.5 uL of ultra pure water (188.75 uL was added for this experiment accidentally)
 - ii. 20 uL of 5x phusion master mix
 - iii. 2 uL 10 mM dNTP
 - iv. 1.25 uL of primer A
 - v. 1.25 uL of primer B
 - vi. 2 uL of DNA template (from step 1)
 - d. Chill mixture for 15 minutes
 - e. Add 1 uL Phusion DNA polymerase with chilled pipette tip
 - f. Transfer 20 uL of sample to each PCR tube with chilled pipette tip
6. PCR cycle
 - a. 95 C for 6 minutes
 - b. 98 C for 30 seconds
 - c. 98 C for 10 seconds
 - d. 52 C to 68 C for 30 seconds
 - i. Sample 1: column 2 @ 52.4C
 - ii. Sample 2: column 4 @ 54.6C
 - iii. Sample 3: column 5 @ 56.5C

- iv. Sample 4: column 7 @ 61.5C
- v. Sample 5: column 9 @ 65.5C
- e. 72 C for 1:45 seconds (30 sec per kb)
- f. Goto step b four times
- g. 98 C for 10 seconds
- h. 67 C for 30 seconds
- i. 72 C for 1:45 seconds
- j. Goto step f twenty nine times
- k. 72 C for 10 minutes
- l. 4C forever
- m. End