

Colony PCR flu operon

8/5/2010

Modifications from 8/4/2010 are highlighted in red

Prepare DNA Template

1. Take a stab of the frozen stock of *E. coli* K12 from the IGEM -80C frozen stock (not too much because too much template will lead to unspecific replication) and resuspend it in 250 uL of DI water (Unspecific amplification so template concentration is too high)
 - a. Add 200 uL of ultra pure water to DNA template already prepared
2. Mix primers to a 50 mM concentration
 - a. Multiply the nmol concentration from the IDT printout by 20
 - b. This is the volume of resuspension in ultrapure water in uL
 - i. Forward flu primer: $24.5 \text{ nmol} \times 20 = 490 \text{ uL}$ ultra pure water
 - ii. Reverse flu primer $46.9 \text{ nmol} \times 20 = 938 \text{ uL}$ ultra pure water
 - c. Done Previously
3. Turn on PCR machine to instant incubate at 98 C for the first initial denaturation step
4. Place PCR tubes for run on ice
5. PCR reaction KEEP EVERYTHING ON ICE
 - a. Vortex all tubes before starting to make sure everything is well mixed
 - b. For a single reaction mixture (50 uL)
 - i. 36.275 uL Ultra pure water
 - ii. 10 uL of 5x phusion master mix
 - iii. 1 uL 10 mM dNTP
 - iv. 0.625 uL of primer A (0.5 uM concentration from 40 uM stock)
 - v. 0.625 uL of primer B (0.5 uM concentration from 40 uM stock)
 - vi. 1 uL of DNA template (from step 1)
 - vii. 0.5 uL Phusion DNA polymerase
 - c. Combine all ingredients but polymerase for 3 samples in master mix (150 ul)
 - i. 108.825 uL of ultra pure water
 - ii. 30 uL of 5x phusion master mix
 - iii. 3 uL 10 mM dNTP
 - iv. 1.875 uL of primer A
 - v. 1.875 uL of primer B
 - vi. 3 uL of DNA template (from step 1)
 - d. Chill mixture for 15 minutes
 - e. Add 1.5 uL Phusion DNA polymerase with chilled pipette tip
 - f. Transfer 50 uL of sample to each PCR tube with chilled pipette tip
6. PCR cycle
 - a. 98 C for 30 seconds
 - b. 98 C for 10 seconds
 - c. 65 C to 68 C for 30 seconds
 - i. Sample 6: column 6 @ 64.8C
 - ii. Sample 7: column 7 @ 66C
 - iii. Sample 8: column 9 @ 67.4C

- d. 72 C for 1:45 seconds (30 sec per kb)
- e. Goto step b four times
- f. 98 C for 10 seconds
- g. 67 C for 30 seconds
- h. 72 C for 1:45 seconds
- i. Goto step f twenty nine times
- j. 72 C for 10 minutes
- k. 4C forever
- l. End