

# Biofilm Formation Experiment

## Background

Naphthenic acids (NA's) are the major contaminant in oil sands production. It is reported that naphtheinic acids (NA's) are best degraded in a biofilm. We want to enhance the biofilm formation ability of bacteria that are capable of degrading NA's in tailing waters. To give the bacteria the ability to form biofilms we plan on introducing a plasmid with antigen 43 that will make the cells "frizzy".

For this experiment we are testing the native ability of bacteria to form biofilms (without the antigen 43 plasmid). Three strains will be tested for biofilm formation: *Pseudomonas putida* (oil sands), *Pseudomonas fluorescens* (oil sands), *Pseudomonas putida* KT2440 and *E. coli* K12. The two *Pseudomonas* strains labeled oilsands were given to us and shown to grow in oil sands tailing water. We will use these *Pseudomonas* strains as model organisms that can degrade NAs. These *Pseudomonas* strains are known to form biofilms and we want to characterize the biofilm formation ability. The *E. coli* K12 and *Pseudomonas putida* KT2440 strain will be used as reference samples.

We will be modeling the adhesion assay from the following paper:

A rapid screening procedure to identify mini-Tn10 insertion mutants of *Escherichia coli* K-12 with altered adhesion properties. FEMS microbiology letters [0378-1097] Genevaux yr:1996 vol:142 iss:1 pg:27

and we will also be adapting our protocol from Alex's lab procedure

## Materials

### Strains

- *E. coli* K12
- *P. putida* (oil sands)
- *P. fluorescens* (oil sands)
- *P. putida* KT2440

### Containers

- 96 well sterile microplate
- 18-15 mL falcon tubes

### Media

- Approx. 12 mL LB media
  - 8 mL for overnight growth
  - 4 mL for microplate inoculation
- Approx. 10 mL bushnell-haas media with glucose carbon source
  - 6mL for overnight growth

- 4 mL for microplate inoculation
- Approx. 3 mL of M9 minimal media with glucose
  - 2 mL for overnight growth
  - 1 mL for microplate inoculation
- Approx. 10 mL bushnell-haas media with glucose and casamino acids
  - 6 mL for overnight growth
  - 4 mL for microplate inoculation
- Approx 3 mL of M9 minimal media with glucose and casamino acids
  - 2 mL for overnight growth
  - 1 mL for microplate inoculation
- Approx. 10 mL bushnell-haas media with cyclohexane carboxylic acid (CHCA)
  - 6 mL for overnight growth
  - 4 mL for microplate inoculation
- Approx. 10 mL bushnell-haas media with cyclohexane carboxylic acid and pH adjusted to 9
  - 6 mL for overnight growth
  - 4 mL for microplate inoculation
- Approx. 13 mL of 0.1% crystal violet staining solution
- Approx 13 mL of 200 proof ethanol
- 35 mL PBS buffer for washing plate

\*\*\*For media recipes see iGEM 2010 media recipes section of the notebook

## ***Procedure***

1. Start an overnight culture of the following strains in 2 mL of of specified media and grow up at 30C in the shaker

Strain	Media
E. coli K12	LB
	M9 + glucose
	M9 + glucose +casamino acids
P. putida KT2440	LB
	BH +glucose
	BH+glucose+casamino acids
	BH+CHCA
	BH+CHCA pH9
P. putida oilsands	LB
	BH +glucose
	BH+glucose+casamino acids
	BH+CHCA
	BH+CHCA pH9
P. fluorescens oilsands	LB
	BH +glucose
	BH+glucose+casamino acids
	BH+CHCA
	BH+CHCA pH9

2. After cultures have been allowed to grow up for specified amount of time above, measure the OD<sub>600</sub> of each culture in the microplate reader
3. Using the OD<sub>600</sub> make a concentrated cell suspension for the microplate inoculation according to the following formula
  - a. Volume for resuspension (μL) = 500 x OD<sub>600</sub>*
  - b. This should give a cell suspension w/ OD600 ~2, which will be used to inoculate the microplate.*
  - c. Following this procedure ensures that each well will be inoculated with the same number of cells*

To make cell suspension

1. Take 1 mL of culture
  2. Spin down at 12000 rpm for 2 minutes
  3. Pipet off supernatant
  4. Resuspend in calculated amount of media from above formula
4. Add 200 uL of media and 2 uL of cell suspension to each microplate well according to the layout below (a polystyrene microplate should be used)

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	DI H2O	LB	LB	M9 glu	M9 glu	M9 glu CAA	M9 glu CAA	BH CHCA	BH CHCA	BH CHCA pH 9	BH CHCA pH 9	DI H2O
<b>B</b>	DI H2O	LB E. coli K12	LB P. fluorescen s oilsands	M9 glucose E. coli K12	BH glucose P. fluorescen s oilsands	M9 Glu. CAA E. coli K12	BH Glu. CAA P. fluorescen s oilsands	BH CHCA P. putida KT244 0	BH CHCA P. putida oilsands	BH CHCA pH 9 P. fluorescens oilsands	DI H2O	DI H2O
<b>C</b>	DI H2O	LB E. coli K12	LB P. fluorescen s oilsands	M9 glucose E. coli K12	BH glucose P. fluorescen s oilsands	M9 Glu. CAA E. coli K12	BH Glu. CAA P. fluorescen s oilsands	BH CHCA P. putida KT244 0	BH CHCA P. putida oilsands	BH CHCA pH 9 P. fluorescens oilsands	DI H2O	DI H2O
<b>D</b>	DI H2O	LB E. coli K12	LB P. fluorescen s oilsands	M9 glucose E. coli K12	BH glucose P. fluorescen s oilsands	M9 Glu. CAA E. coli K12	BH Glu. CAA P. fluorescen s oilsands	BH CHCA P. putida KT244 0	BH CHCA P. putida oilsands	BH CHCA pH 9 P. fluorescens oilsands	DI H2O	DI H2O
<b>E</b>	DI H2O	LB P. putida KT24 40	LB P. putida oilsands	BH glucose putida KT244 0	BH glucose P. putida oilsands	BH Glu. CAA P. putida KT244 0	BH Glu. CAA P. putida oilsands	BH CHCA P. fluoresc ens oilsands	BH CHCA pH 9 P. putida KT2440	BH CHCA pH 9 P. putida oilsands	DI H2O	DI H2O
<b>F</b>	DI H2O	LB P. putida KT24 40	LB P. putida oilsands	BH glucose putida KT244 0	BH glucose P. putida oilsands	BH Glu. CAA P. putida KT244 0	BH Glu. CAA P. putida oilsands	BH CHCA P. fluoresc ens oilsands	BH CHCA pH 9 P. putida KT2440	BH CHCA pH 9 P. putida oilsands	DI H2O	DI H2O
<b>G</b>	DI H2O	LB P. putida KT24 40	LB P. putida oilsands	BH glucose putida KT244 0	BH glucose P. putida oilsands	BH Glu. CAA P. putida KT244 0	BH Glu. CAA P. putida oilsands	BH CHCA P. fluoresc ens oilsands	BH CHCA pH 9 P. putida KT2440	BH CHCA pH 9 P. putida oilsands	DI H2O	DI H2O
<b>H</b>	DI H2O	LB	LB	BH glu	BH glu	BH glu CAA	BH glu CAA	BH CHCA	BH CHCA	BH CHCA pH 9	BH CHCA pH 9	DI H2O

Key:

M9=M9 minimal media

BH=Bushnell-Haas minimal media

CAA=Casamino acids

CHCA=cyclohexane carboxylic acid

- Allow the plate to incubate (WITHOUT SHAKING OR ANYTYPE OF MOTION) for 48 hours at 30C
- Measure the OD600 on the plate reader after the 48 hour period
  - This step will be used to normalize the faster growing cultures against the slower growing ones
- Pipet out liquid culture with multiwell pipette

- a. Include the blank wells for all of the CV staining steps
8. Wash plate with 200  $\mu$ L of PBS
9. Stain cells with 200  $\mu$ L of 0.1% crystal violet stain for 15 minute
10. Wash cells two more times with PBS
11. Let the plates dry completely
12. Add 200  $\mu$ L of ethanol to each well, cover plate and let it sit for 10 minute
13. Briefly mix the contents of each well
14. Measure the adsorption at OD600 to quantify the biofilm formation

## ***Results***

Comparing the adsorption readings of the OD600 of the cultures after biofilm growth (OD) vs the OD 600 after the crystal violet staining CV-OD600 should give us the biofilm formation comparison between the three strains. A plot of CV-OD600/OD600 vs the different strains will be made.