

Colony PCR flu operon
8/11/2010

Prepare DNA Template

1. Mix primers to a 50 mM concentration
 - a. Multiply the nmol concentration from the IDT printout by 20
 - b. This is the volume of resuspension in ultrapure water in uL
 - i. Forward flu primer: $24.5 \text{ nmol} \times 20 = 490 \text{ uL}$ ultra pure water
 - ii. Reverse flu primer $46.9 \text{ nmol} \times 20 = 938 \text{ uL}$ ultra pure water
 - c. Done Previously
2. Turn on PCR machine to instant incubate at 98 C for the first initial denaturation step
3. Place PCR tubes for run on ice
4. PCR reaction KEEP EVERYTHING ON ICE
 - a. Vortex all tubes before starting to make sure everything is well mixed
 - b. For a single reaction mixture (50 uL)
 - i. 36.275 uL Ultra pure water
 - ii. 10 uL of 5x phusion master mix
 - iii. 1 uL 10 mM dNTP
 - iv. 0.625 uL of primer A (0.5 uM concentration from 40 uM stock)
 - v. 0.625 uL of primer B (0.5 uM concentration from 40 uM stock)
 - vi. 1 uL of DNA template (flu operon gel purified product)
 - vii. 0.5 uL Phusion DNA polymerase
 - c. Combine all ingredients but polymerase for 3 samples in master mix (150 ul)
 - i. 108.825 uL of ultra pure water
 - ii. 30 uL of 5x phusion master mix
 - iii. 3 uL 10 mM dNTP
 - iv. 1.875 uL of primer A
 - v. 1.875 uL of primer B
 - vi. 3 uL of DNA template (flu operon gel purified product)
 - d. Chill mixture for 15 minutes
 - e. Add 1.5 uL Phusion DNA polymerase with chilled pipette tip
 - f. Transfer 50 uL of sample to each PCR tube with chilled pipette tip
5. PCR cycle
 - a. 98 C for 30 seconds
 - b. 98 C for 10 seconds
 - c. 52 to 68 C for 30 seconds
 - i. Sample 16 in column 8 at 63.8 C
 - ii. Sample 17 in column 9 at 65.5 C
 - iii. Sample 18 in column 12 at 68 C
 - d. 72 C for 1:45 seconds (30 sec per kb)
 - e. Goto step b 29 times
 - f. 72 C for 10 minutes
 - g. 4C forever
 - h. End