

## **Enhanced Arginase production: rocF**

### **Purpose and Justification:**

*Bacillus subtilis* produces urease, which catalyses the hydrolysis of urea into ammonium and carbonate. Since the cell wall of the bacteria is negatively charged, the bacteria draw cations from the environment, including  $\text{Ca}^{2+}$ , to deposit on their cell surface. The  $\text{Ca}^{2+}$  ions subsequently react with the  $\text{CO}_3^{2-}$  ions, leading to the precipitation of  $\text{CaCO}_3$  at the cell surface.

In order for *B. subtilis* to fill up cracks in concrete, enhanced production of calcium carbonate must be achieved.

Previous experiments involving up-regulating *ureA*, *ureB* and *ureC* in *B. subtilis* have not led to an increase in urease production. This could be due to yet unidentified genes that are involved in the process. Therefore we are looking at another strategy, of increasing urea hydrolysis indirectly by increasing arginine and arginase production. Arginase breaks down arginine to urea and ornithine, leading to an increase of urea inside the cell. We believe that in turn the urea itself will increase urease production. By targeting this pathway, it will hopefully be able to activate additional unknown pathways and enzymes that are involved in calcium carbonate production.

We plan to produce two BioBricks, *SR1*, which will enhance arginine production (please see the *SR1* BioBrick cloning strategy), and *rocF*, which will enhance arginase production.

The naturally occurring gene *rocF* codes for the protein arginase.

## Modelling:

We have done computational modelling of our system. The graph below shows that the final output of the pathway, carbonate, increases as time progresses.

Details of this model can be found on the wiki.




Fig 1. A graph showing the final output of the biochemical pathway, carbonate, increases over time

## Construction:

### By Gibson method

Parts:

1. 01bbprefix - The standard BioBrick prefix.
2. 03pspac\_oid – The Pspacoid promoter including the oid operator, repressed by *lacI*. Allows us to test the brick by adding IPTG. This sequence is from part Bba\_K174004 on the Parts Registry.
3. 05rocF\_cds - The coding sequence for *rocF*. This coding sequence is from GenBank ([X81802.1](#)).
4. 06dbTerm09 - A double terminator. This is part Bba\_B0014 on the parts registry.
5. 07bb\_suffix – The standard BioBrick suffix.




Fig 2. Screenshot of the *rocF* BioBrick components in Sequencer

The BioBrick is going to be constructed using the Gibson method (*Enzymatic assembly of DNA molecules up to several hundred kilobases*, Gibson et al.). Overlapping PCR products will be joined together to form the final construct shown in Fig 2.

The *rocF* coding sequence has two EcoR1 restriction sites which must be removed if the construct is to be BioBrick compatible (see Fig 3.). In order to do this, we will PCR the *rocF* coding sequence as three overlapping products and join them using the Gibson method. Our primers will have single base changes which preserve the amino acid sequence but make the coding sequence BioBrick-compatible (see Fig 4. , 5. and 6.).




Fig 3. A restriction map of the *rocF* coding sequence, showing two EcoR1 sites we must remove

Sequencher - (000001-000001) (including TBS)

File Edit Select Contig Sequence View Window Help

Overview Cut Map Ruler Find

Residue: 227 to 227

```

1      TCITAAATGT TGAGGTGGAA TCACAGATGG ATAAAACGAT
      L N V E V E S Q M D K T I
41     TTCGGTTATT GGAATGCCAA TGGATTTAGG ACAAGCACGA
      S V I G M P M D L G Q A R
81     CGCGGAGTGG ATATGGGCCG GAGTGCCATC CGGTACGCTC
      R G V D M G P S A I R Y A H
121    ATCTGATCGA GAGGCTGTCA GACATGGGGT ATACGGTTGA
      L I E R L S D M G Y T V E
161    AGATCTCGGT GACATTCCGA TCAATCCGGA AAAAAACAAA
      D L G D I P I N R E K I K
201    AATGACGAGG AACTGAAAAA CCTGAATTCG GTTTTGGCGG
      N D E E L K N L N S V L A G
241    GAAATGAAAA ACTCGGCCAA AAGGTCAACA AAGTCATTGA
      N E K L A Q K V N K V I E
281    AGAGAAAAAA TTCGCCCTTG TCCTGGGCGG TGACCACAGT
      E K K F P L V L G G D H S
321    ATTGCGATCG GCACGGTTGC AGGCACAGCG AAGCATTACG
      I A I G T L A G T A K H Y D
361    ATAATCTCGG CGTCATCTGG TATGACGCGC ACGGGCGATT
      N L G V I W Y D A H G D L
401    GAATACACTT GAAACTTCAC CATCGGGCAA TATTCACGGC
      N T L E T S P S G N I H G
441    ATGCCGCTCG CGGTACGCCT AGGCATTGGC CACGAGTCAC
      M P L A V S L G I G H E S L
481    TGGTTAACCT TGAAGGCTAC GCGCCTAAAA TCAAACCGGA
      V N L E G Y A P K I K P E
521    AAACGTCGTC ATCATTGGCG CCCGGTCACT TGATGAAGGG
      N V V I I G A R S L D E G
561    GAGCGCAAGT ACATTAAGGA AAGCGGCATG AAGGTGTACA
      E R K Y I K E S G M K V Y T
601    CAATGCACGA AATCGATCGT CTGGCATGA CAAAGTTCAT
      M H E I D R L G M T K V I
641    TGAAGAAACC CTTGATTATT TATCAGCATG TGATGGCGTC
      E E T L D Y L S A C D G V
681    CATCTGAGCC TTGATCTGGA CGGACTTGAT CCGAACGAGC
  
```

Sequencher - (000001-000001)

File Edit Select Contig Sequence View Window Help

Overview Cut Map Ruler Find Show Experiments

Residue: 227 to 227

```

1      TCTTAAATGT TGAGGTGGAA TCACAGATGG ATAAAACGAT
      L N V E V E S Q M D K T I
41     TTCGGTTATT GGAATGCCAA TGGATTTAGG ACAGCACGA
      S V I G M P M D L G Q A R
81     CGCGGAGTGG ATATGGGCCG GAGTGCCATC CGGTACGCTC
      R G V D M G P S A I R Y A H
121    ATCTGATCGA GAGGCTGTCA GACATGGGGT ATACGGTTGA
      L I E R L S D M G Y T V E
161    AGATCTCGGT GACATTCCGA TCAATCCGGA AAAAAACAAA
      D L G D I P I N R E K I K
201    AATGACGAGG AACTGAAAAA CCTGAATTCG GTTTTGGCGG
      N D E E L K N L N S V L A G
241    GAAATGAAAA ACTCGGCCAA AAGGTCAACA AAGTCATTGA
      N E K L A Q K V N K V I E
281    AGAGAAAAAA TTCGCCCTTG TCCTGGGCGG TGACCACAGT
      E K K F P L V L G G D H S
321    ATTGCGATCG GCACGGTTGC AGGCACAGCG AAGCATTACG
      I A I G T L A G T A K H Y D
361    ATAATCTCGG CGTCATCTGG TATGACGCGC ACGGGCGATT
      N L G V I W Y D A H G D L
401    GAATACACTT GAAACTTCAC CATCGGGCAA TATTCACGGC
      N T L E T S P S G N I H G
441    ATGCCGCTCG CGGTACGCCT AGGCATTGGC CACGAGTCAC
      M P L A V S L G I G H E S L
481    TGGTTAACCT TGAAGGCTAC GCGCCTAAAA TCAAACCGGA
      V N L E G Y A P K I K P E
521    AAACGTCGTC ATCATTGGCG CCCGGTCACT TGATGAAGGG
      N V V I I G A R S L D E G
561    GAGCGCAAGT ACATTAAGGA AAGCGGCATG AAGGTGTACA
      E R K Y I K E S G M K V Y T
601    CAATGCACGA AATCGATCGT CTGGCATGA CAAAGTTCAT
      M H E I D R L G M T K V I
641    TGAAGAAACC CTTGATTATT TATCAGCATG TGATGGCGTC
      E E T L D Y L S A C D G V
681    CATCTGAGCC TTGATCTGGA CGGACTTGAT CCGAACGAGC
  
```

Fig 4. Removal of the first EcoR1 restriction site by a single base change





Fig 5. Removal of the second EcoRI restriction site by a single base change




Fig 6. A restriction map of the *rocF* coding sequence with the two single base changes

Below (Fig 7.), the primers we are using are shown. In total there are 6 pairs:

1. Primer 1-P1 and Primer 2-P1 – Promoter PCR product. Primer 1-P1 includes standard BioBrick prefix on 5' end. Product overlaps with product 2, below.
2. Primer 1-S1 and Primer 2-S1 – First fragment of the *rocF* RBS and coding sequence. This product overlaps with product 1 above, and product 3, below.
3. Primer 3-S2 and Primer 4-S2 – Second fragment of the *rocF* coding sequence. This product overlaps with product 2 above, and product 4, below.
4. Primer 5-S3 and Primer 6-S3 – Third and final fragment of the *rocF* coding sequence. This product overlaps with product 3 above, and product 5, below.
5. Primer 1-T1 and Primer 2-T1 – Double terminator PCR product. Primer 2-T1 includes standard BioBrick suffix on 5' end. Product overlaps with product 4, above.
6. Primer 1-V1 and Primer 2-V1 – Vector primers. Used for cloning into a BioBrick-compatible plasmid, again using the Gibson method (see *Procedure*, below). Primer 1-V1 is the standard BioBrick suffix, and Primer 2-V1 is the standard BioBrick prefix.




Fig 7. The primers used to PCR the overlapping promoter, *rocF* coding sequence and RBS, and double terminator products






Fig 8. Removal of the first EcoRI restriction site in the product by single base change in Primer 2-S1 and Primer 3-S2





Fig 9. Removal of the second EcoR1 restriction site in the product by single base change in Primer 4-S2 and Primer 5-S3

### Map of the *rocF* Biobrick



## Procedure:

1. Take 6 PCR tubes and label them 1-6.
2. The part which would be amplified in each tube are as following:

| Tube | Part to be amplified                        |
|------|---|
| 1    | BioBrick-compatible Plasmid vector          |
| 2    | Pspacoid Promoter                           |
| 3    | 1 <sup>st</sup> fragment of <i>rocF</i> CDS |
| 4    | 2 <sup>nd</sup> fragment of <i>rocF</i> CDS |
| 5    | 3 <sup>rd</sup> fragment of <i>rocF</i> CDS |
| 6    | Double terminator                           |

3. Add 1 µl of the plasmid (template DNA) extracted from *E. coli* DH5α into all the tubes.
4. Add the PCR reagents as mentioned below in the tubes:

| Number       | Reagents                   | Volume  |
|--------------|----------------------------|---------|
| 1            | Distilled H <sub>2</sub> O | 27.5 µl |
| 2            | 5x PCR Buffer              | 10 µl   |
| 3            | Nucleotide dNTPs           | 1 µl    |
| 4            | Forward primer             | 5 µl    |
| 5            | Backward primer            | 5 µl    |
| 6            | Phusion Polymerase         | 0.5 µl  |
| 7            | Template DNA               | 1 µl    |
| Total Volume | -                          | 50 µl   |

Now, the forward and reverse which should go into each PCR tube is:

| Tube | Part to be amplified                        | Forward Primer | Reverse Primer |
|------|---|----------------|----------------|
| 1    | Plasmid vector                              | Primer 1- V1   | Primer 2-V1    |
| 2    | Pspacoid Promoter                           | Primer 1- P1   | Primer 2- P1   |
| 3    | 1 <sup>st</sup> fragment of <i>rocF</i> CDS | Primer 1- S1   | Primer 2- S1   |
| 4    | 2 <sup>nd</sup> fragment of <i>rocF</i> CDS | Primer 3- S2   | Primer 4- S2   |
| 5    | 3 <sup>rd</sup> fragment of <i>rocF</i> CDS | Primer 5- S3   | Primer 6- S3   |
| 6    | Double terminator                           | Primer 1- T1   | Primer 2- T1   |

5. Set the conditions as mentioned below for the thermocycler as following after putting the PCR tubes into the PCR wells present on the thermocycler:

| Step             | Condition                      | Time         |
|------------------|--------------------------------|--------------|
| Initialise       | 98°C                           | 30 seconds   |
| Denature*        | 98°C                           | 10 seconds   |
| Anneal*          | Depends on the size of the DNA | 20 seconds   |
| Extension*       | 72°C                           | 30 seconds   |
| Extension finish | 72°C                           | 5-10 minutes |
| Hold             | 4°C                            | N/A          |

\* Remember to set these steps for 30 cycles.

6. After the completion of the PCR reaction, run Gel electrophoresis (0.7 % agarose) to check whether the PCR reaction has worked.

7. Extract the band (Refer to the gel extraction protocol) of the right size after referring to the DNA ladder in the first lane.

| Tube | Part to be amplified                        | Size of the fragment |
|------|---|----------------------|
| 1    | Plasmid vector                              | 2072 bp              |
| 2    | Pspacoid Promoter                           | +106 bp              |
| 3    | 1 <sup>st</sup> fragment of <i>rocF</i> CDS | + 246 bp             |
| 4    | 2 <sup>nd</sup> fragment of <i>rocF</i> CDS | +597 bp              |
| 5    | 3 <sup>rd</sup> fragment of <i>rocF</i> CDS | +125 bp              |
| 6    | Double terminator                           | +116 bp              |

8. Set up a single PCR tube for the final step for the Gibson cloning. This tube would contain the 6 amplified DNA fragments from above and would be annealed in the subsequent steps.

9. For the Gibson cloning method, the following volumes of the chemicals are required:

(A) 1.33X Master Mix

| Reagents                 | Volumes |
|--------------------------|---------|
| 5X isothermal buffer     | 100 µl  |
| T5 exonuclease 1.0 U/ µl | 2 µl    |
| Phusion DNA pol 2 U/ µl  | 6.25 µl |
| Taq DNA ligase 40 U/ µl  | 50 µl   |

|                  |           |
|------------------|-----------|
| H <sub>2</sub> O | 216.75 µl |
| Total Volume*    | 375 µl    |

\*This volume makes 25 aliquots of 15 µl each.

(B) 5X Isothermal Reaction Buffer

| Reagents                | Volumes |
|-------------------------|---------|
| 25% PEG-8000            | 0.75 g  |
| 500 mM Tris-HCl pH 7.5  | 1.5 ml  |
| 50 mM MgCl <sub>2</sub> | 75 µl   |
| 50 mM DTT               | 150 µl  |
| 1 mM dATP               | 30 µl   |
| 1 mM dTTP               | 30 µl   |
| 1 mM dCTP               | 30 µl   |
| 1 mM dGTP               | 30 µl   |
| 5 mM NAD                | 300 µl  |
| H <sub>2</sub> O        | 105 µl  |
| Total Volume            | 3 ml    |

(C) Final Volume in the PCR Tube

| Reagents  | Volumes |
|---|---------|
| Template DNA (consisting of all 6 parts which were amplified from the PCR reaction mentioned above) | 5 µl    |
| 1.33 X Master Mix   | 15 µl   |
| Total Volume  | 20 µl   |

\* Remember to do this step very quickly as everything has to be done on ice.

10. Incubate the cells at 50°C for 60 minutes for the reaction to take place.

## Integration:

The fragment is to be cloned into the BioBrick-compatible version of pGFPrrnB and integrated into the *Bacillus subtilis* 168 chromosome at *amyE*.

## Testing and Characterisation:

### Selection for integration

To select for integration of the plasmid into the chromosome, *B. subtilis* will be tested for the ability to hydrolyse starch. Integration of the BioBrick will be done by homologous recombination at the *amyE* gene, therefore destroying endogenous expression of amylase. Colonies that are not able to break down starch on agar plate will be selected and cultured for further test. Colonies that do not contain the integrated BioBrick will be able to hydrolyse starch, therefore forming a white halo around the colony as iodine interacts with starch to form blue colour.

### Characterisation of *rocF*

The aim of this experiment is to determine whether the *rocF* BioBrick increases arginase production.

### Materials Required

- Plate consisting of *Bacillus subtilis* 168 colonies.
- Flame (streaking) Loop
- LB media consisting arginine and ampicillin
- Auto pipette
- Bursen Burner
- Universal Tube

### Procedure

- Perform the experiment using aseptic technique.
- Transfer *B. subtilis* 168 colonies into universal tubes containing 5 ml of LB media and allowed to grow overnight at 37° C.
- Transfer 1 ml of the overnight culture to another universal tube containing 4 ml of the following media:
  1. Control (1) - LB media
  2. Control (2) - LB media with 10 mM of arginine
  3. Control (3) - LB media plus *B. subtilis* 168
  4. Test (1) - LB media with 10 mM of arginine plus *B. subtilis* 168
  5. Test (2) - LB media with 10 mM of arginine plus *B. subtilis* 168
- Incubate the culture at 37° C with shaking.
- Record the pH at every 30 min interval. Use 20 ul of the culture and measure the pH.

## **Expected results**

Arginase breaks down arginine to ornithine. Increasing production of arginase in the cells will cause a faster breakdown of arginine, thus faster changes in pH level.

1. Control (1) - No change in pH
2. Control (2) - No change in pH
3. Control (3) - Increase in pH, however will be lower than test 1 and test 2.
4. Test (1) - Increase in pH
5. Test (2) - Increase in pH