

Competent cells



Swiss Federal Institute of Technology Lausanne iGEM team 2010, Dupont Thibault, Gerweck Nadia, Grädel Nadine, Helfer Jonas, Lisowski Wiktor, Monnot Gwennaëlle, Perrudet Christian, Richter Solange, Varicchio Stefano, Vokali Efthymia, Déneraud Nicolas, Gubelmann Carine, Niederholtmeyer Henrike Marie, Opota Onya, Deplancke Bart, Maerkl Sebastian.

Introduction

This protocol provides information on how to make *Asaia* competent cells for approximately 10 transformations. This process takes 3 days.

Remarks: be careful with the competent cells. Thaw them only when you need them and always keep the cells on ice. Do not shake or pipette competent cells too vigorously. We rec-

ommend to use microtubes for each transformation, so you don't have to refreeze the competent *Asaia* after each use.

MATERIAL

To make competent cells, you need :

- ◆ *Asaia* liquid culture
- ◆ An optical density reader
- ◆ A -80°C freezer
- ◆ A centrifuge
- ◆ Liquid nitrogen
- ◆ 10% cold Glycerol
- ◆ 1mM HEPES at pH5

PROTOCOL

Day 1

1. Pick-up a colony from a plate and make an overnight pre-culture in 2 ml of Glycerol medium (GLY medium).

Day 2

2. Transfer 1 ml of the pre-culture in 49 ml of GLY medium in a 500 ml flask. Put it at 30°C overnight.

Day 3

3. Dilute liquid culture with the ratio 1:11 in GLY medium, e.g. 20 ml of culture with 200 ml of GLY medium.
4. Incubate with aeration until cells reach early log

phase (optical density at 550 nm between 0.5 and 0.8).

5. Transfer culture into a 15 ml centrifuge tube.
6. Incubate them on ice for 15 minutes. *After this point, it's very important to keep the cells cold!*
7. Sediment them at $2'700\text{g}$ for 10 minutes at 4°C .
8. Remove the supernatant.
9. Re-suspend the pellet with 10ml of 1mM HEPES at pH5.
10. Pellet the bacteria at $2'700\text{g}$ for 10 min at 4°C .
11. Remove the supernatant.
12. Redo step 8 to 10 once
13. Re-suspend the cells in 5 ml of cold 10% Glycerol.
14. Sediment them at $2'700\text{g}$ for 10 minutes at 4°C .
15. Remove the supernatant.
16. Re-suspend the cells in 0.65 ml of cold 10% Glycerol.
17. Fill microtubes with 65 μl of competent cells.
18. Snap freeze the tubes in liquid nitrogen to freeze them.
19. Put all tubes in the -80°C fridge.